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Explorer

Are instruments used to examine and measure. They should never be employed to remove dental tartar. This category of instruments can basically be divided into two different families: periodontal and explorer.

Periodontal explorer

Are essentially used to measure the depth of the periodontal groove and pockets. Their tips are rounded to prevent damage to tissues. These instruments are equipped with coded reference notches that enable precise measurements to be taken. The instrument is inserted under the marginal gingiva. It is held parallel to the surface of the tooth, and along its main axis. Measurement is taken by gently pushing the tip downwards until it reaches the groove or pocket. A soft, elastic resistance (somewhat like a taut strip of rubber) is felt when the bottom is reached. At this point, the depth can be measured on the notches in the instrument. A hard, solid resistance during measurement indicates the presence of tartar.

Explorer probes

Are used to detect the presence of tartar, caries and surface irregularities on the tooth, protrusions and open junctions in reconstructions, decalcified areas, and lesions due to caries. Basically, there are various types of explorer probes:

- a) **The curved explorer** (Fig. 03), which is used to detect tartar and caries. Excellent for revealing tartar in deep pockets or bifurcations. The thin tip provides good tactile sensitivity. It also enables the instrument to adapt properly to the sub-gingival environment and to any kind of tooth surface. The back of the tip can be placed against the base of the pocket without lacerating the epithelial junction.
- b) **"Orban"-type explorer** (fig. 17), which is principally used to reveal subgingival tartar in tight, deep pockets since it adapts to these areas without causing inconvenient dislocation of the tissues. The back of the tip can be placed against the bottom of the pocket without lacerating the epithelial junction.
- c) **Explorer with shepherd's hook** (fig. 23), which is principally used to search for caries. Generally, it is too rigid and thick to provide good tactile sensitivity. Its shape prevents it from adapting to interproximal areas, bifurcations and deep pockets. Not recommended for revealing tartar by direct application of the tip on the epithelium because the junction could be perforated.

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Root canal instruments

Pluggers - These instruments are used in endodontics and to measure the depth of the root canal for subsequent devitalisation when executing root canal fillings.

Spreaders - The heated tip is used to melt gutta-percha.

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All rights reserved. Reproduction of the images and texts in this catalog by any means, whole or in part, without Den Tag permission is prohibited.